

Antibodies Product Datasheet

Product Name:	P5CR1 Antibody
Catalog Number:	C51165
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal P5CR1
Concentration:	1 mg/ml
Gene Name:	PYCR1
Protein Name:	Pyrroline-5-carboxylate reductase 1, mitochondrial (P5C reductase 1) (P5CR 1) (EC 1.5.1.2)
Other Name:	Pyrroline-5-carboxylate reductase 1, mitochondrial (P5C reductase 1) (P5CR 1) (EC 1.5.1.2)
Human Swiss Prot No.:	P32322
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of P5CR1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year