
Antibodies Product Datasheet

Product Name:	Retinoic Acid Receptor α/β (Phospho-Ser96) Antibody
Catalog Number:	A8222
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Retinoic Acid Receptor α/β (Phospho-Ser96)
Concentration:	1 mg/ml
Gene Name:	RARA NR1B1
Protein Name:	Retinoic Acid Receptor α/β (Phospho-Ser96)
Other Name:	Retinoic acid receptor alpha (RAR-alpha) (Nuclear receptor subfamily 1 group B member 1)
Human Swiss Prot No.:	P10276
Formulation:	Liquid in PBS containing 50% glycerol, and 0.207% sodium azide.
Specificity:	This antibody detects endogenous levels of Retinoic Acid Receptor α/β (Phospho-Ser96) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year