

Antibodies Product Datasheet

Product Name:	KEPI Antibody
Catalog Number:	C11380
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP1R14C. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	PPP1R14C
Protein Name:	Protein phosphatase 1 regulatory subunit 14C
Other Name:	PPP1R14C; KEPI; Protein phosphatase 1 regulatory subunit 14C; Kinase-enhanced PP1 inhibitor; PKC-potentiated PP1 inhibitory protein; Serologically defined breast cancer antigen NY-BR-81
Human Swiss Prot No.:	Q8TAE6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KEPI Polyclonal Antibody detects endogenous levels of KEPI protein.
Storage and Stability:	-20°C/1 year