

Antibodies Product Datasheet

Product Name:	CRYZL1 Antibody
Catalog Number:	C50133
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CRYZL1
Concentration:	1 mg/ml
Gene Name:	CRYZL1 4P11
Protein Name:	Quinone oxidoreductase-like protein 1 (EC 1.-.-) (Protein 4P11) (Quinone oxidoreductase homolog 1) (QOH-1) (Zeta-crystallin homolog)
Other Name:	
Human Swiss Prot No.:	O95825
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CRYZL1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year