
Antibodies Product Datasheet

Product Name:	TAL-1 Antibody
Catalog Number:	R12-3586
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	TAL-1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAL-1.
Concentration:	1 mg/ml
Gene Name:	TAL1
Other Name:	SCL; Stem cell protein; T-cell acute lymphocytic leukemia 1 (SCL); T-cell acute lymphocytic leukemia-1 protein; T-cell leukemia/lymphoma-5 protein; TAL1; TCL5
Human Swiss Prot No.:	P17542
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TAL-1 Antibodyantibody detects endogenous levels of TAL-1.
Storage and Stability:	Stable at -20°C for at least 1 year.