

Antibodies Product Datasheet

Product Name:	Rb (Phospho-Ser249) Antibody
Catalog Number:	A0982
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Retinoblastoma around the phosphorylation site of Ser249. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	RB1
Protein Name:	Retinoblastoma-associated protein
Other Name:	RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110
Human Swiss Prot No.:	P06400
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Rb (S249) Polyclonal Antibody detects endogenous levels of Rb protein only when phosphorylated at S249.
Storage and Stability:	-20°C/1 year