

Antibodies Product Datasheet

Product Name:	Ferritin Light Chain mouse mAb
Catalog Number:	N0987
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant full length of human ferritin light chain protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ftl
Protein Name:	
Other Name:	Ferritin L chain ; Ferritin L subunit ; Ferritin light chain ; Ferritin light polypeptide ; Ferritin light polypeptide like 3 ; FRIL ; FRIL_HUMAN ; FTL ; L apoferritin ; LFTD ; MGC71996 ; NBIA 3 ; NBIA3.
Human Swiss Prot No.:	P02792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects recombinant ferritin proteins.
Storage and Stability:	-20°C/1 year