

Antibodies Product Datasheet

Product Name:	CREB-2 Antibody
Catalog Number:	B7017
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATF4. AA range:212-261
Concentration:	1 mg/ml
Gene Name:	ATF4
Protein Name:	Cyclic AMP-dependent transcription factor ATF-4
Other Name:	ATF4; CREB2; TXREB; Cyclic AMP-dependent transcription factor ATF-4; cAMP-dependent transcription factor ATF-4; Activating transcription factor 4; Cyclic AMP-responsive element-binding protein 2; CREB-2; cAMP-responsive element-binding prot
Human Swiss Prot No.:	P18848
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CREB-2 Polyclonal Antibody detects endogenous levels of CREB-2 protein.
Storage and Stability:	-20°C/1 year