Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Skp1 p19 Antibody

Catalog Number: C0323

Reactivity: Human, Mouse, Rat, Monkey

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000. Not yet tested in other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human SKP1A/p19. AA range:41-90

Concentration: 1 mg/ml

Gene Name: SKP1

Protein Name: S-phase kinase-associated protein 1

Other Name: SKP1; EMC19; OCP2; SKP1A; TCEB1L; S-phase kinase-associated

protein 1; Cyclin-A/CDK2-associated protein p19; Organ of Corti protein

2; OCP-2; Organ of Corti protein II; OCP-II; RNA polymerase II

elongation factor-like protein; SIII; Transcr

Human Swiss Prot No.: P63208

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Skp1 p19 Polyclonal Antibody detects endogenous levels of Skp1 p19

protein.

Storage and Stability: -20°C/1 year