

Antibodies Product Datasheet

Product Name:	Caspase 2 (p13, Cleaved-Gly334) Antibody
Catalog Number:	L0479
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Caspase 2 (p13, Cleaved-Gly334)
Concentration:	1 mg/ml
Gene Name:	CASP2 ICH1 NEDD2
Protein Name:	Caspase2
Other Name:	Caspase-2 (CASP-2) (EC 3.4.22.55) (Neural precursor cell expressed developmentally down-regulated protein 2) (NEDD-2) (Protease ICH-1) [Cleaved into: Caspase-2 subunit p18; Caspase-2 subunit p13; Caspase-2 subunit p12]
Human Swiss Prot No.:	P42575
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Caspase 2 (p13, Cleaved-Gly334, protein was cleaved amino acid sequence between 333-334)
Storage and Stability:	-20°C/1 year