
Antibodies Product Datasheet

Product Name:	Collagen XI α 1 Antibody
Catalog Number:	C12214sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Collagen XI α 1 Polyclonal
Concentration:	1 mg/ml
Gene Name:	COL11A1 COL6
Protein Name:	Collagen alpha-1(XI) chain
Other Name:	Collagen alpha-1(XI) chain
Human Swiss Prot No.:	P12107
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Collagen XI α 1.
Storage and Stability:	-20°C/1 year