
Antibodies Product Datasheet

Product Name:	p53 Monoclonal Antibody(6C4)
Catalog Number:	N1316
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:2000 IHC: 1:200 IF: 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of p53
Concentration:	
Gene Name:	TP53
Protein Name:	Cellular tumor antigen p53
Other Name:	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Human Swiss Prot No.:	P04637
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous p53 proteins.
Storage and Stability:	-20°C/1 year