
Antibodies Product Datasheet

| | |
|---------------------------------|---|
| Product Name: | Anti-Calponin-1 [ABT-CALP] |
| Catalog Number: | V0013 |
| Clone ID: | ABT-CALP |
| Reactivity: | Human |
| Applications: | IHC |
| Application Dilutions | IHC 1:100-200 |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| UniProt: | P51911 |
| Gene Name: | CNN1 |
| Other Name: | Calponin H1, smooth muscle |
| Isotype | Mouse IgG |
| Localization: | N/A |
| Positive Tissue Control: | N/A |
| Immunogen | Synthetic peptide |
| Purification | Immunogen affinity purified |
| Formulation: | PBS, pH7.4, 50% glycerol, 0.03%Proclin300 |
| Specificity: | Anti-Calponin-1 [ABT-CALP] detects endogenous levels of total Calponin-1 protein. |
| Storage and Stability: | -20°C/1 year |