

Antibodies Product Datasheet

Product Name:	WWP1 Antibody
Catalog Number:	C51198
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal WWP1
Concentration:	1 mg/ml
Gene Name:	WWP1
Protein Name:	NEDD4-like E3 ubiquitin-protein ligase WWP1 (EC 6.3.2.-) (Atrophin-1-interacting protein 5) (AIP5) (TGIF-interacting ubiquitin ligase 1) (Tiul1) (WW domain-containing protein 1)
Other Name:	NEDD4-like E3 ubiquitin-protein ligase WWP1 (EC 6.3.2.-) (Atrophin-1-interacting protein 5) (AIP5) (TGIF-interacting ubiquitin ligase 1) (Tiul1) (WW domain-containing protein 1)
Human Swiss Prot No.:	Q9H0M0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of WWP1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year