

Antibodies Product Datasheet

Product Name:	RASL2 Antibody
Catalog Number:	C45800
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	RASA4 CAPRI GAPL KIAA0538
Protein Name:	Ras GTPase-activating protein 4 (Calcium-promoted Ras inactivator) (Ras p21 protein activator 4) (RasGAP-activating-like protein 2)
Other Name:	
Human Swiss Prot No.:	O43374
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	RASL2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year