

Antibodies Product Datasheet

Product Name:	DNA-PKCS Antibody
Catalog Number:	B0909
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DNA-PK. AA range:2022-2071
Concentration:	1 mg/ml
Gene Name:	PRKDC
Protein Name:	DNA-dependent protein kinase catalytic subunit
Other Name:	PRKDC; HYRC; HYRC1; DNA-dependent protein kinase catalytic subunit; DNA-PK catalytic subunit; DNA-PKcs; DNPK1; p460
Human Swiss Prot No.:	P78527
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DNA-PKCS Polyclonal Antibody detects endogenous levels of DNA-PKCS protein.
Storage and Stability:	-20°C/1 year