

Antibodies Product Datasheet

Product Name:	Adiponectin Antibody
Catalog Number:	L0131
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Adiponectin . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	ADIPOQ
Protein Name:	Adiponectin
Other Name:	ADIPOQ; ACDC; ACRP30; APM1; GBP28; Adiponectin; 30 kDa adipocyte complement-related protein; Adipocyte complement-related 30 kDa protein; ACRP30; Adipocyte; C1q and collagen domain-containing protein; Adipose most abundant gene transcript 1
Human Swiss Prot No.:	Q15848
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Adiponectin Polyclonal Antibody detects endogenous levels of Adiponectin protein.
Storage and Stability:	-20°C/1 year