

Antibodies Product Datasheet

Product Name:	hnRNP C1/2 Antibody
Catalog Number:	B1041
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human hnRNP C1/2. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	HNRNPC
Protein Name:	Heterogeneous nuclear ribonucleoproteins C1/C2
Other Name:	HNRNPC; HNRPC; Heterogeneous nuclear ribonucleoproteins C1/C2; hnRNP C1/C2
Human Swiss Prot No.:	P07910
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	hnRNP C1/2 Polyclonal Antibody detects endogenous levels of hnRNP C1/2 protein.
Storage and Stability:	-20°C/1 year