

Antibodies Product Datasheet

Product Name:	Anti-Human Claudin 3 (ABT509) Mouse Monoclonal Antibody
Catalog Number:	V0262
Clone ID:	ABT509
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CLDN3
Gene Name:	O15551
Other Name:	C7orf1, Claudin-3, Claudin3, CLD3_HUMAN, CLDN 3, Cldn3, Clostridium perfringens enterotoxin receptor 2, CPE R2, CPE receptor 2, CPE-R 2, CPE-receptor 2, CPETR 2, CPETR2, HRVP 1, HRVP1, Rat ventral prostate 1 like protein, Rat ventral prostate.1 protein homolog, RVP1, Ventral prostate.1 like protein, Ventral prostate.1 protein homolog
Isotype	IgG2b, kappa
Localization:	Colon/ Appendix
Positive Tissue Control:	Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Claudin 3 protein.
Storage and Stability:	-20°C/1 year