
Antibodies Product Datasheet

Product Name:	TTL11 Antibody
Catalog Number:	C46918
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TTL11
Concentration:	1 mg/ml
Gene Name:	TLL11 C9orf20
Protein Name:	Tubulin polyglutamylase TLL11 (EC 6.3.2.1) (Tubulin--tyrosine ligase-like protein 11)
Other Name:	
Human Swiss Prot No.:	Q8NHH1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human TTL11
Storage and Stability:	-20°C/1 year