

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | SAP30BP Antibody  |
| <b>Catalog Number:</b>        | C50052  |
| <b>Reactivity:</b>            | Human, Mouse  |
| <b>Applications:</b>          | WB  |
| <b>Dilutions:</b>             | WB 1:500-2000   |
| <b>Modification:</b>          | Unmodified/Total  |
| <b>Source:</b>                | Polyclonal, Rabbit, IgG   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human SAP30BP  |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | SAP30BP HCNGP HTRG HTRP   |
| <b>Protein Name:</b>          | SAP30-binding protein (Transcriptional regulator protein HCNGP)   |
| <b>Other Name:</b>            |   |
| <b>Human Swiss Prot No.:</b>  | Q9UHR5  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of SAP30BP at Human, Mouse  |
| <b>Storage and Stability:</b> | Store at -20° C/1 year  |