

Antibodies Product Datasheet

Product Name:	BM28 Monoclonal Antibody
Catalog Number:	N0069
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human BM28 expressed in E. Coli.
Concentration:	
Gene Name:	MCM2
Protein Name:	DNA replication licensing factor MCM2
Other Name:	MCM2; BM28; CCNL1; CDCL1; KIAA0030; DNA replication licensing factor MCM2; Minichromosome maintenance protein 2 homolog; Nuclear protein BM28
Human Swiss Prot No.:	P49736
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	BM28 Monoclonal Antibody detects endogenous levels of BM28 protein.
Storage and Stability:	-20°C/1 year