

Antibodies Product Datasheet

Product Name:	CTR2 Antibody
Catalog Number:	C51169
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-ternal CTR2
Concentration:	1 mg/ml
Gene Name:	SLC7A2 ATRC2 CAT2
Protein Name:	Low affinity cationic amino acid transporter 2 (CAT-2) (CAT2) (Solute carrier family 7 member 2)
Other Name:	Low affinity cationic amino acid transporter 2 (CAT-2) (CAT2) (Solute carrier family 7 member 2)
Human Swiss Prot No.:	P52569
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTR2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year