

Antibodies Product Datasheet

Product Name:	Uba5 Antibody
Catalog Number:	C19417
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human UBA5. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	UBA5
Protein Name:	Ubiquitin-like modifier-activating enzyme 5
Other Name:	UBA5; UBE1DC1; Ubiquitin-like modifier-activating enzyme 5; Ubiquitin-activating enzyme 5; ThiFP1; UFM1-activating enzyme; Ubiquitin-activating enzyme E1 domain-containing protein 1
Human Swiss Prot No.:	Q9GZZ9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Uba5 Polyclonal Antibody detects endogenous levels of Uba5 protein.
Storage and Stability:	-20°C/1 year