

Antibodies Product Datasheet

Product Name:	SAMHD1 (Phospho-Thr592) Antibody
Catalog Number:	A9491
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human SAMHD1 (Thr592)
Concentration:	1 mg/ml
Gene Name:	SAMHD1 MOP5
Protein Name:	SAMHD1 (Thr592)
Other Name:	SAM domain and HD domain-containing protein 1 (EC 3.1.4.-) (Dendritic cell-derived IFNG-induced protein) (DCIP) (Monocyte protein 5) (MOP-5)
Human Swiss Prot No.:	Q9Y3Z3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat SAMHD1 (phospho-Thr592)
Storage and Stability:	-20°C/1 year