

Antibodies Product Datasheet

| | |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name: | Cleaved-MMP-15 (Y132) Antibody |
| Catalog Number: | L0323 |
| Reactivity: | Human, Monkey |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Cleaved Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human MMP15. AA range:113-162 |
| Concentration: | 1 mg/ml |
| Gene Name: | MMP15 |
| Protein Name: | Matrix metalloproteinase-15 |
| Other Name: | MMP15; Matrix metalloproteinase-15; MMP-15; Membrane-type matrix metalloproteinase 2; MT-MMP 2; MTMMP2; Membrane-type-2 matrix metalloproteinase; MT2-MMP; MT2MMP; SMCP-2 |
| Human Swiss Prot No.: | P51511 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Cleaved-MMP-15 (Y132) Polyclonal Antibody detects endogenous levels of fragment of activated MMP-15 protein resulting from cleavage adjacent to Y132. |
| Storage and Stability: | -20°C/1 year |