

Antibodies Product Datasheet

Product Name:	CaV α 2 δ 3 Antibody
Catalog Number:	C40053
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CaV α 2 δ 3
Concentration:	1 mg/ml
Gene Name:	CACNA2D3
Protein Name:	Voltage-dependent calcium channel subunit alpha-2/delta-3 (Voltage-gated calcium channel subunit alpha-2/delta-3) [Cleaved into: Voltage-dependent calcium channel subunit alpha-2-3; Voltage-dependent calcium channel subunit delta-3]
Other Name:	Voltage-dependent calcium channel subunit alpha-2/delta-3 (Voltage-gated calcium channel subunit alpha-2/delta-3) [Cleaved into: Voltage-dependent calcium channel subunit alpha-2-3; Voltage-dependent calcium channel subunit delta-3]; YM3539
Human Swiss Prot No.:	Q8IZS8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CaV α 2 δ 3 protein(A210) detects endogenous levels of CaV α 2 δ 3
Storage and Stability:	-20°C/1 year