
Antibodies Product Datasheet

Product Name:	MYPT1 (Phospho-Ser507) Antibody
Catalog Number:	A9049
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human MYPT1 (Ser507)
Concentration:	1 mg/ml
Gene Name:	PPP1R12A MBS MYPT1
Protein Name:	MYPT1 (Ser507)
Other Name:	Protein phosphatase 1 regulatory subunit 12A (Myosin phosphatase-targeting subunit 1) (Myosin phosphatase target subunit 1) (Protein phosphatase myosin-binding subunit)
Human Swiss Prot No.:	O14974
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat MYPT1 (phospho-Ser507)
Storage and Stability:	-20°C/1 year