

Antibodies Product Datasheet

Product Name:	TRPV1 Antibody
Catalog Number:	C46874
Reactivity:	Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 440-520
Concentration:	1 mg/ml
Gene Name:	TRPV1 VR1
Protein Name:	Transient receptor potential cation channel subfamily V member 1 (TrpV1) (Capsaicin receptor) (Osm-9-like TRP channel 1) (OTRPC1) (Vanilloid receptor 1)
Other Name:	
Human Swiss Prot No.:	Q8NER1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	TRPV1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year