
Antibodies Product Datasheet

Product Name:	TTR mouse Monoclonal Antibody(1D7)
Catalog Number:	N1668
Reactivity:	Human
Applications:	IF/ICC,WB,IHC-p
Dilutions:	IF: 1:50-200 WB 1:500-2000,IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of TTR
Concentration:	1 mg/ml
Gene Name:	TTR PALB
Protein Name:	Transthyretin (ATTR) (Prealbumin) (TBPA)
Other Name:	Transthyretin (ATTR) (Prealbumin) (TBPA)
Human Swiss Prot No.:	P02766
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous TTR protein
Storage and Stability:	-20°C/1 year