

Antibodies Product Datasheet

Product Name:	GBAS Antibody
Catalog Number:	C50795
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GBAS
Concentration:	1 mg/ml
Gene Name:	GBAS NIPSNAP2
Protein Name:	Protein NipSnap homolog 2 (NipSnap2) (Glioblastoma-amplified sequence)
Other Name:	
Human Swiss Prot No.:	O75323
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GBAS at Human
Storage and Stability:	Store at -20° C/1 year