
Antibodies Product Datasheet

Product Name:	p57 Antibody
Catalog Number:	C10659
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDKN1C. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	CDKN1C
Protein Name:	Cyclin-dependent kinase inhibitor 1C
Other Name:	CDKN1C; KIP2; Cyclin-dependent kinase inhibitor 1C; Cyclin-dependent kinase inhibitor p57; p57Kip2
Human Swiss Prot No.:	P49918
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p57 Polyclonal Antibody detects endogenous levels of p57 protein.
Storage and Stability:	-20°C/1 year