

Antibodies Product Datasheet

Product Name:	PLEKHF1 Antibody
Catalog Number:	C49386
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PLEKHF1
Concentration:	1 mg/ml
Gene Name:	PLEKHF1 APPD LAPF ZFYVE15
Protein Name:	Pleckstrin homology domain-containing family F member 1 (PH domain-containing family F member 1) (Lysosome-associated apoptosis-inducing protein containing PH and FYVE domains) (Apoptosis-inducing protein) (PH and FYVE domain-containing protein 1) (Phafin
Other Name:	
Human Swiss Prot No.:	Q96S99
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PLEKHF1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year