

Antibodies Product Datasheet

Product Name:	YTDC1 Antibody
Catalog Number:	C51193
Reactivity:	Human, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal YTDC1
Concentration:	1 mg/ml
Gene Name:	YTHDC1 KIAA1966 YT521
Protein Name:	YTH domain-containing protein 1 (Putative splicing factor YT521)
Other Name:	YTH domain-containing protein 1 (Putative splicing factor YT521)
Human Swiss Prot No.:	Q96MU7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of YTDC1 at Human,Rat
Storage and Stability:	Store at -20° C/1 year