

Antibodies Product Datasheet

Product Name:	RNF15 Antibody
Catalog Number:	C10031
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRIM38. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	TRIM38
Protein Name:	Tripartite motif-containing protein 38
Other Name:	TRIM38; RNF15; RORET; Tripartite motif-containing protein 38; RING finger protein 15; Zinc finger protein RoRet
Human Swiss Prot No.:	O00635
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RNF15 Polyclonal Antibody detects endogenous levels of RNF15 protein.
Storage and Stability:	-20°C/1 year