

## Antibodies Product Datasheet

<b>Product Name:</b>	OAS2 Antibody
<b>Catalog Number:</b>	C30992
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from the N-terminal region of human OAS2. AA range:61-110
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	OAS2
<b>Protein Name:</b>	2'-5'-oligoadenylate synthase 2
<b>Other Name:</b>	OAS2; 2'-5'-oligoadenylate synthase 2; (2-5')oligo(A) synthase 2; 2-5A synthase 2; p69 OAS / p71 OAS; p69OAS / p71OAS
<b>Human Swiss Prot No.:</b>	P29728
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	OAS2 Polyclonal Antibody detects endogenous levels of OAS2 protein.
<b>Storage and Stability:</b>	-20°C/1 year