

Antibodies Product Datasheet

Product Name:	Cytokeratin 18 (PT1876) mouse mAb
Catalog Number:	N1836
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgM, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cytokeratin 18
Concentration:	
Gene Name:	KRT18 CYK18 PIG46
Protein Name:	Keratin, type I cytoskeletal 18 (Cell proliferation-inducing gene 46 protein) (Cytokeratin-18) (CK-18) (Keratin-18) (K18)
Other Name:	
Human Swiss Prot No.:	P05783
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Cytokeratin 18. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year