
Antibodies Product Datasheet

Product Name:	N-CoR Antibody
Catalog Number:	C0358
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NCoR1. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	NCOR1
Protein Name:	Nuclear receptor corepressor 1
Other Name:	NCOR1; KIAA1047; Nuclear receptor corepressor 1; N-CoR; N-CoR1
Human Swiss Prot No.:	O75376
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	N-CoR Polyclonal Antibody detects endogenous levels of N-CoR protein.
Storage and Stability:	-20°C/1 year