

## Antibodies Product Datasheet

<b>Product Name:</b>	LIN28A (Phospho-Ser200) Antibody
<b>Catalog Number:</b>	A9736
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho peptide around human LIN28A (Ser200)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	LIN28A CSDD1 LIN28 ZCCHC1
<b>Protein Name:</b>	LIN28A (Ser200)
<b>Other Name:</b>	Protein lin-28 homolog A (Lin-28A) (Zinc finger CCHC domain-containing protein 1)
<b>Human Swiss Prot No.:</b>	Q9H9Z2
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human LIN28A (phospho-Ser200)
<b>Storage and Stability:</b>	-20°C/1 year