

Antibodies Product Datasheet

Product Name:	TCP-1 η Antibody
Catalog Number:	C19085
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCT7. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	CCT7
Protein Name:	T-complex protein 1 subunit eta
Other Name:	CCT7; CCTH; NIP7-1; T-complex protein 1 subunit eta; TCP-1-eta; CCT-eta; HIV-1 Nef-interacting protein
Human Swiss Prot No.:	Q99832
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TCP-1 η Polyclonal Antibody detects endogenous levels of TCP-1 η protein.
Storage and Stability:	-20°C/1 year