

## Antibodies Product Datasheet

<b>Product Name:</b>	FoxO4 (Acetyl-Lys189) Antibody
<b>Catalog Number:</b>	D12076
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000, ELISA 1:10000-20000
<b>Modification:</b>	Acetyl Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthetic Acetyl peptide from human protein at AA range: 189
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	FOXO4 AFX AFX1 MLLT7
<b>Protein Name:</b>	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
<b>Other Name:</b>	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
<b>Human Swiss Prot No.:</b>	P98177
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	The antibody detects endogenous FoxO4 when Acetyl occurs at Lys189
<b>Storage and Stability:</b>	-20°C/1 year