
Antibodies Product Datasheet

Product Name:	Anti-Caldesmon [ABT-CALD]
Catalog Number:	V0012
Clone ID:	ABT-CALD
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	Q05682
Gene Name:	CALD1
Other Name:	CAD; CALD1; CDM; L-caldesmon; Non-muscle caldesmon
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Caldesmon [ABT-CALD] detects endogenous levels of total Caldesmon protein.
Storage and Stability:	-20°C/1 year