

Antibodies Product Datasheet

Product Name:	ZNF397 Antibody
Catalog Number:	C11371
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF397. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	ZNF397
Protein Name:	Zinc finger protein 397
Other Name:	ZNF397; ZNF47; ZSCAN15; Zinc finger protein 397; Zinc finger and SCAN domain-containing protein 15; Zinc finger protein 47
Human Swiss Prot No.:	Q8NF99
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ZNF397 Polyclonal Antibody detects endogenous levels of ZNF397 protein.
Storage and Stability:	-20°C/1 year