

Antibodies Product Datasheet

Product Name:	ELOVL3 Antibody
Catalog Number:	C15619
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELOVL3. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	ELOVL3
Protein Name:	Elongation of very long chain fatty acids protein 3
Other Name:	ELOVL3; CIG30; Elongation of very long chain fatty acids protein 3; 3-keto acyl-CoA synthase ELOVL3; Cold-inducible glycoprotein of 30 kDa; ELOVL fatty acid elongase 3; ELOVL FA elongase 3
Human Swiss Prot No.:	Q9HB03
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ELOVL3 Polyclonal Antibody detects endogenous levels of ELOVL3 protein.
Storage and Stability:	-20°C/1 year