

Antibodies Product Datasheet

Product Name:	Cytokeratin 8 (Phospho-Ser74) Antibody
Catalog Number:	P12-1017
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Cytokeratin 8 (Phospho-Ser74) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Cytokeratin 8 around the phosphorylation site of Ser74.
Concentration:	1 mg/ml
Gene Name:	KRT8
Other Name:	Keratin, type II cytoskeletal 8; Keratin-8; Type-II keratin Kb8; CK-8
Human Swiss Prot No.:	P05787
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Cytokeratin 8 (Phospho-Ser74) antibody detects endogenous levels of Cytokeratin 8 only when phosphorylated at Ser74.
Storage and Stability:	Stable at -20°C for at least 1 year.