

Antibodies Product Datasheet

Product Name:	Dyrk1A Antibody
Catalog Number:	C11002
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DYR1A. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	DYRK1A
Protein Name:	Dual specificity tyrosine-phosphorylation-regulated kinase 1A
Other Name:	DYRK1A; DYRK; MNB; MNBH; Dual specificity tyrosine-phosphorylation-regulated kinase 1A; Dual specificity YAK1-related kinase; HP86; Protein kinase minibrain homolog; MNBH; hMNB
Human Swiss Prot No.:	Q13627
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dyrk1A Polyclonal Antibody detects endogenous levels of Dyrk1A protein.
Storage and Stability:	-20°C/1 year