

Antibodies Product Datasheet

Product Name:	p164-RhoGEF Antibody
Catalog Number:	C11524
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARHGEF17. AA range:431-480
Concentration:	1 mg/ml
Gene Name:	ARHGEF17
Protein Name:	Rho guanine nucleotide exchange factor 17
Other Name:	ARHGEF17; KIAA0337; TEM4; Rho guanine nucleotide exchange factor 17; 164 kDa Rho-specific guanine-nucleotide exchange factor; p164-RhoGEF; p164RhoGEF; Tumor endothelial marker 4
Human Swiss Prot No.:	Q96PE2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p164-RhoGEF Polyclonal Antibody detects endogenous levels of p164-RhoGEF protein.
Storage and Stability:	-20°C/1 year