

Antibodies Product Datasheet

Product Name:	Ephrin-A4 Antibody
Catalog Number:	C0181
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EFNA4. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	EFNA4
Protein Name:	Ephrin-A4
Other Name:	EFNA4; EPLG4; LERK4; Ephrin-A4; EPH-related receptor tyrosine kinase ligand 4; LERK-4
Human Swiss Prot No.:	P52798
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ephrin-A4 Polyclonal Antibody detects endogenous levels of Ephrin-A4 protein.
Storage and Stability:	-20°C/1 year