
Antibodies Product Datasheet

Product Name:	Vitronectin Antibody
Catalog Number:	C30305
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human VTN. AA range:50-100
Concentration:	1 mg/ml
Gene Name:	VTN
Protein Name:	Vitronectin
Other Name:	VTN; Vitronectin; VN; S-protein; Serum-spreading factor; V75
Human Swiss Prot No.:	P04004
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Vitronectin Polyclonal Antibody detects endogenous levels of Vitronectin protein.
Storage and Stability:	-20°C/1 year