

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | ADCK5 Antibody |
| Catalog Number: | C11156 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human ADCK5. AA range:151-200 |
| Concentration: | 1 mg/ml |
| Gene Name: | ADCK5 |
| Protein Name: | Uncharacterized aarF domain-containing protein kinase 5 |
| Other Name: | ADCK5; Uncharacterized aarF domain-containing protein kinase 5 |
| Human Swiss Prot No.: | Q3MIX3 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | ADCK5 Polyclonal Antibody detects endogenous levels of ADCK5 protein. |
| Storage and Stability: | -20°C/1 year |